

Audiology

How to use your personal listener



Your personal listener can be used when watching TV or having a conversation with someone.



Insert a 9 volt battery in the back of the personal listener.



Connect the headphones where you see this symbol on the side.



Put the headphones on and turn the dial until it is loud enough.

Any questions?

If you have any questions about the information contained in this leaflet, then please contact the Audiology department to help you.

If you would like this information in other formats, such as...

- Another language
- Large print
- Audio

then please contact a member of the Audiology department.

Audiology Department Box 94

Addenbrooke's Hospital, Cambridge University Hospitals NHS Foundation Trust, Hills Road, Cambridge, CB2 2QQ

Tel: 01223 217797

Fax: 01223 586912

Text phone: 01223 274494

Document History

Authors	Audiology Department
Department	Cambridge University Hospitals NHS Foundation Trust, Hills Road, Cambridge, CB2 0QQ www.cuh.org.uk
Contact Number	01223 596360
Publish/ Review Date	August 2010 / August 2013
File Name	How_to_use_your_personal_listener
Version number/ Ref	1/PIN2652